Subject: Fwd: Filter Wheel

Date: Tuesday, December 21, 2021 at 11:44:57 AM Eastern Standard Time

From: Stephen Hope
To: Massimo Robberto

Attachments: 12162_5.png, PastedGraphic-5.tiff

External Email - Use Caution

Regards,

PastedGraphic-5.tiff

Stephen C. Hope

Electronic Engineering/Software Development



Johns Hopkins University

Instrument Development Group Bloomberg Center for Physics & Astronomy 3701 San Martin Drive Baltimore, MD 21218

Telephone: +1(410) 516 4161 Email: <u>stephen.hope@idg.jhu.edu</u>

Web: http://idg.jhu.edu

Begin forwarded message:

From: stephen.hope@idg.jhu.edu
Subject: Fwd: Filter Wheel

Date: October 14, 2021 at 11:40:47 AM EDT **To:** Jonathan Hoover < 12jonboy12@gmail.com >

Regards,

PastedGraphic-5.tiff

Stephen C. Hope Electronic Engineering/Software Development



Johns Hopkins University Instrument Development Group Bloomberg Center for Physics & Astronomy 3701 San Martin Drive Baltimore, MD 21218

Telephone: +1(410) 516 4161 Email: <u>stephen.hope@idg.jhu.edu</u>

Web: http://idg.jhu.edu

Begin forwarded message:

From: stephen.hope@idg.jhu.edu

Subject: Filter Wheel

Date: May 21, 2021 at 12:03:14 PM EDT

To: Jonathan Hoover < 12jonboy12@gmail.com >

Here are filter wheel commands to send to the PCM.... The **stored procedures** only need to be sent once in the lifetime of the drive... You should have the ability to send these, but will not have to unless we replaced a drive.

The **operating commands** allow you to choose a particular filter by asking the drive to execute a stored procedure. This should be integrated into your code. You also ought to query the wheels to make sure they are where you think they are....

Note: all procedures include the command to route the drive packet via the PCM so it should be really simple to integrate.

BTW the camera is no facing the filter wheel.

<u>FILTER WHEEL STORED PROCEDURES</u> - These need to be sent to the motor controllers the first time they are used... i.e if you replace a motor controller you should reinitialize with these commands...

WHEEL 1

Initial drive settings (velocity, steps ratio, current), plus homing routine — this runs automatically on power up.

~@,9600_8N1T2000,/1s0m23l23h0j32V4000v2500n0P0P100z0M500e1R

Position A1

~@,9600_8N1T2000,/1s1A46667R

Position A2

~@,9600_8N1T2000,/1s2A62222R

Position A3

~@,9600_8N1T2000,/1s3A77778R

Position A4

~@,9600_8N1T2000,/1s4A0R

Position A5

~@,9600_8N1T2000,/1s5A15555R

Position A6

~@,9600_8N1T2000,/1s6A31111R

WHEEL 2

```
Initial drive settings (velocity, steps ratio, current), plus homing routine — this runs automatically on power up.
```

~@,9600_8N1T2000,/2s0m23l23h0j32V4000v2500n0P0P100z0M500e1R

Position B1

~@,9600_8N1T2000,/2s1A46667R

Position B2

~@,9600_8N1T2000,/2s2A62222R

Position B3

~@,9600_8N1T2000,/2s3A77778R

Position B4

~@,9600_8N1T2000,/2s4A0R

Position B5

~@,9600_8N1T2000,/2s5A15555R

Position B6

~@,9600_8N1T2000,/2s6A31111R

OPERATING COMMANDS:

Homing occurs automatically on power up.... Or it can be forced using:

```
~@,9600 8N1T2000,/1e0R (wheel 1)
```

~@,9600_8N1T2000,/2e0R (wheel 2)

Positions 1 through six are commanded as follows:

```
~@,9600_8N1T2000,/1e1R (Position A1)
```

~@,9600_8N1T2000,/2e1R (Position B1)

~@,9600_8N1T2000,/1e2R (Position A2)

~@,9600_8N1T2000,/2e2R (Position B2)

~@,9600_8N1T2000,/1e3R (Position A3)

~@,9600_8N1T2000,/2e3R (Position B3)

~@,9600_8N1T2000,/1e4R (Position A4)

~@,9600_8N1T2000,/2e4R (Position B4)

```
~@,9600_8N1T2000,/1e5R (Position A5)
~@,9600_8N1T2000,/2e5R (Position B5)

~@,9600_8N1T2000,/1e6R (Position A6)
~@,9600_8N1T2000,/2e6R (Position B6)

QUERY COMMANDS:

~@,9600_8N1T2000,1/?0 (returns wheel 1, current step count)
~@,9600_8N1T2000,/2?0 (returns wheel 2, current step count)
```

Example response: /0`15555

When the wheel is in position (i.e. not moving) the step count should be:

- Position 1: 46667
 Position 2: 62222
 Position 3: 77778
 Position 4: 0
 Position 5: 15555
 Position 6: 31111
- ~@,9600_8N1T2000,1/?0 (returns wheel 1, sensor status) ~@,9600_8N1T2000,/2?0 (returns wheel 2, sensor status)

Example Response: /0'14

When the wheel is in position (i.e. not moving) the sensor status should be:

Position 1: 14
Position 2: 14
Position 3: 14
Position 4: 13
Position 5: 14
Position 6: 14

Regards,

PastedGraphic-5.tiff



Johns Hopkins University Instrument Development Group Bloomberg Center for Physics & Astronomy 3701 San Martin Drive Baltimore, MD 21218

Telephone: +1(410) 516 4161 Email: <u>stephen.hope@idg.jhu.edu</u>

Web: http://idg.jhu.edu